

Product Name:	Little Dino Kid's Toothbrush
Materials:	PP, PBT, Silicone
Dimensions:	116 x 25 x 13mm
Colours:	Baby Blue/Pale Pink
Recommended Age:	2 years+

Product Description:

Say goodbye to the struggle of brushing teeth with the fun Little Dino Kid's Toothbrush! Designed to make brushing comfortable for little ones, it features dense, ultra-soft bristles that effectively clean teeth without hurting gums. The rounded brush head is perfectly sized for small mouths, making it easy for babies to use without discomfort. The raised bumps on the brush head gently clean residue on your child's tongue. Featuring a cute, anti-slip and soft dinosaur-shaped handle with a sturdy base to stand stably on the sink, this toothbrush is easy for toddlers to grip, helping them develop those all-important dental hygiene habits and transition to regular toothbrushes. The simple design enables easy cleaning – simply boil to sterilise!

Cleaning and Care:

Wash before each use. We recommend hand washing with warm soapy water, then rinsing and drying thoroughly. Boil in water for 2-3 minutes to sterilise. Do not use bleach-based cleaners or tablets to clean this product. Store product in a cool, dry place and avoid direct sunlight. We recommend replacing your child's toothbrush every 3-4 months.

Instructions:

Wet the bristles with clean water. Apply the desired amount of toddler-safe toothpaste*. Let your child grab the brush to clean their teeth and gums using gentle circular motions. Ensure that all tooth surfaces are cleaned. Rinse and remove all foam from child's mouth. Place the cleaned brush on the sink.

*Follow your dentist's guidelines around the use of toothpaste for children. Adult supervision is required during use, and children should be guided on proper brushing techniques.

Caution:

Check product condition regularly. Replace if this product shows any signs of damage. Avoid storing near sharp objects. This product is not a toy. Do not use this product for anything other than its intended purpose. Avoid allowing children to chew on or pull the toothbrush to prevent bristles from deforming or falling out. If the bristles show signs of wear, splitting, or deformation, replace the toothbrush immediately. Adult supervision is required. Avoid playing with the toothbrush while brushing to prevent accidents. Keep out of reach of children and pets when not in use. NOT suitable for pets.

PLEASE NOTE: Brush head is detachable. If it slips off, simply reattach by pressing back onto the handle.